

POSTER PRESENTATION

Open Access

Isometric exercise in cardiac magnetic resonance imaging: an initial experience using fast imaging

Kristian H Mortensen^{1*}, Jennifer A Steeden², Joseph Panzer¹, Andrew M Taylor¹, Vivek Muthurangu¹

From 2011 SCMR/Euro CMR Joint Scientific Sessions
Nice, France. 3-6 February 2011

Introduction

Isometric exercise is a powerful stimulator of the cardiovascular system. Therefore, exercise can be used to unmask subtle cardiovascular changes that are not evident at rest. MRI is a proven method for comprehensive cardiovascular assessment in the resting state. Unfortunately, lack of fast imaging sequences and an appropriate exercise test has prohibited the use of isometric exercise in MRI.

Purpose

We aimed to devise an isometric exercise test for use in MRI in order to comprehensively assess the hemodynamic response to exercise using fast-imaging MRI.

Methods

A sustained isometric biceps exercise test was devised for the MRI scanner. Ten healthy male volunteers (24 to 44 years) underwent three minutes of sustained isometric exercise using a weight load adjusted to 35% of maximum biceps strength. The hemodynamic response

was elucidated during rest, exercise [1, 2 and 3 minutes] and recovery [1, 3 and 5 minutes]. At each time point we obtained 1) aortic flow - spiral prospectively-gated breath-hold phase contrast CMR, 2) LV volumes - real-time radial k-t sensitivity encoding free-breathing CMR and 3) oscillometric upper arm blood pressure. Using this data vascular resistance and compliance was calculated.


Results

All participants completed the isometric exercise test. All data sets were complete. Sustained isometric exercise in the MRI scanner produced a marked hemodynamic response (Table 1). Both vascular and left ventricular measures were changed by the exercise. Most indices were changed at one minute, while all were changed after two and three minutes (Figure 1). The return towards the resting state at recovery was fast (Table 1). There was no evidence of compensatory adaptive hemodynamic mechanisms during the recovery period in the cardiovascular measures.

Table 1 Mean (SD) cardiovascular measures during rest, isometric exercise and recovery

	Rest	Isometric biceps exercise			Recovery		
		1 min	2 min	3 min	1 min	3 min	5 min
Systolic blood pressure (mmHg)	111 (10)	128 (16) †	140 (21) †	145 (17) †	116 (13) †	113 (12)	112 (11)
Heart rate (beats.min ⁻¹)	64 (9)	80 (8) †	81 (9) †	82 (13) †	65 (8)	64 (9)	63 (9)
Stroke volume (mL.m ⁻²)	54 (8)	51 (9) †	49 (9) †	49 (8) †	55 (8)	56 (9)	54 (8)
Cardiac output (L.min ⁻¹ .m ⁻²)	3.4 (0.5)	4.1 (0.8) †	4.0 (0.8) †	4.0 (0.7) †	3.6 (0.4)	3.6 (0.5)	3.4 (0.4)
Compliance (mmHg.mL ⁻¹ .m ⁻¹)	1.3 (0.3)	1.1 (0.2) †	1.0 (0.2) †	0.9 (0.3) †	1.1 (0.2)	1.1 (0.3)	1.2 (0.2)
Resistance (WU.m ²)	24.6 (3.4)	24.7 (4.6)	28.1 (6.6) †	28.1 (4.9) †	24.2 (3.3)	23.8 (3.2)	25.1 (3.4)
End-diastolic volume (mL.m ⁻²)	82 (12)	86 (16) †	85 (15) †	86 (14) †	82 (13)	82 (13)	82 (14)
Ejection fraction (%)	67 (9)	60 (9) †	58 (9) †	58 (6) †	66 (8)	68 (6)	66 (5)

¹UCL Institute of Child Health, London, UK



Conclusion

Isometric biceps exercise testing is highly feasible using a MRI safe methodology, which allows exercise in the supine position within the scanner. The test evokes a hemodynamic response that corresponds to isometric exercise conducted outside the MRI scanner. The hemodynamic response can be comprehensively and safely assessed at any point of time during isometric exercise using this exercise methodology together with fast-imaging MRI.

Author details

¹UCL Institute of Child Health, London, UK. ²UCL Department of Medical Physics & Bioengineering, London, UK.

Published: 2 February 2011

doi:10.1186/1532-429X-13-S1-P386

Cite this article as: Mortensen *et al.*: Isometric exercise in cardiac magnetic resonance imaging: an initial experience using fast imaging. *Journal of Cardiovascular Magnetic Resonance* 2011 **13**(Suppl 1):P386.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at
www.biomedcentral.com/submit

 **BioMed Central**