

Research Article

Validation of copy number variation analysis for next-generation sequencing diagnostics

Jamie M Ellingford^{1,2},
Christopher Campbell¹,
Stephanie Barton¹,
Sanjeev Bhaskar¹,
Saurabh Gupta³,
Rachel L Taylor^{1,2},
Panagiotis I Sergouniotis¹,
Bradley Horn¹,
Janine A Lamb⁴,
Michel Michaelides^{5,6},
Andrew R Webster^{5,6},
William G Newman^{1,2},
Binay Panda³,
Simon C Ramsden¹,
Graeme CM Black^{1,2}

¹Manchester Centre for Genomic Medicine, Central Manchester University Hospitals NHS Foundation Trust, Manchester Academic Health Sciences Centre, St Mary's Hospital, Manchester M13 9WL, UK.

²Institute of Human Development, University of Manchester, Oxford Road, Manchester, M13 9WL, UK.

³Ganit Labs, Bio-IT Centre, Institute of Bioinformatics and Applied Biotechnology, Bangalore, 560100, India.

⁴Institute of Population Health, University of Manchester, Oxford Road, Manchester, M13 9PT, UK.

⁵Moorfields Eye Hospital NHS Foundation Trust, London, EC1V 2PD, UK

⁶UCL Institute of Ophthalmology, Department of Genetics, London, EC1V 9EL, UK.

Corresponding Author: Jamie M Ellingford,
Jamie.ellingford@postgrad.manchester.ac.uk, +44 161 276 8703.

Conflict of Interest: The authors declare no conflict of interest.

Running Title: CNV detection in targeted NGS diagnostics

Abstract

Although a common cause of disease, copy number variants (CNVs) have not routinely been identified from next-generation sequencing (NGS) data in a clinical context. This study aimed to examine the sensitivity, and specificity of a widely used software package, ExomeDepth, to identify CNVs from targeted NGS datasets. We benchmarked the accuracy of CNV detection using ExomeDepth v1.1.16 applied to targeted NGS datasets, through comparison to CNV events detected through whole genome sequencing (WGS) for 25 individuals, and determined the sensitivity and specificity of ExomeDepth applied to these targeted NGS datasets to be 100% and 99.8%, respectively. To define quality assurance metrics for CNV surveillance through ExomeDepth, we undertook simulation of single exon ($n=1000$) and multiple-exon heterozygous deletion events ($n=1749$), determining a sensitivity of 97% ($n=2749$). We identified that the extent of sequencing coverage, the inter- and intra-sample variability in the depth of sequencing coverage, and the composition of analysis regions are all important determinants of successful CNV surveillance through ExomeDepth. We then applied these quality assurance metrics during CNV surveillance for 140 individuals across 12 distinct clinical areas, encompassing over 500 potential rare disease diagnoses. All 140 individuals lacked molecular diagnoses after routine clinical NGS testing, and through application of ExomeDepth we identified 17 CNVs contributing to the cause of a Mendelian disorder. Our findings support the integration of CNV detection using ExomeDepth v1.1.16 with routine targeted NGS diagnostic services for Mendelian disorders. Implementation of this strategy increases diagnostic yields and enhances clinical care.

Key words: Copy Number Variation, Targeted Next-Generation Sequencing, Mendelian disorders, Medical Genetics

Introduction

Molecular diagnostic services available for patients with genetically heterogeneous Mendelian disease have been transformed by the adoption of next-generation DNA sequencing (*NGS*) within the clinical setting.^{1, 2} At present, diagnostic services facilitated by *NGS* are frequently limited to targeted capture techniques, including custom gene panels^{3,4} and whole exome sequencing (*WES*).^{5, 6} These techniques have demonstrated tremendous power to identify rare and private single nucleotide variation and small insertions/deletions underpinning disease onset.

The identification of large structural variants and copy number variants (*CNVs*) encapsulating the regions targeted by *WES* and custom gene panel assays have proved challenging in a clinical context. While whole genome sequencing (*WGS*) techniques have the potential to address this gap in diagnostic *NGS* services,^{7,8} the cost and data burdens remain substantial. Consequently, the application of *CNV* detection algorithms in targeted *NGS* diagnostic services can facilitate immediate improvement in clinical care for individuals with heterogeneous Mendelian disorders. However, such techniques require formal assessment to demonstrate accuracy, reliability and repeatability.

Here, we assess a framework for the implementation of *CNV* detection with targeted *NGS* diagnostic services applied across a range of highly heterogeneous Mendelian disorders.

Methods

Study Design

High coverage targeted NGS data was generated in a United Kingdom Accredited Clinical Laboratory. We applied a CNV detection algorithm to validate the sensitivity for (i) known CNV events, and (ii) simulated CNV events (Figure 1). We assessed a number of factors to determine whether they influenced successful CNV surveillance. We selected two key factors identified from assessments of simulated and known CNVs (inter-sample variability and insufficient coverage) as quality assurance metrics during prospective CNV detection for individuals without molecular diagnoses through clinical NGS testing (Figure 1).

Our analyses included individuals referred for diagnostic testing for four highly heterogeneous disorders where targeted gene panel NGS is a routine diagnostic service, specifically: inherited retinal dystrophies (IRD), congenital cataracts, cardiac disorders and metabolic disorders.

Sequencing & Variant Analysis

Whole Genome and Targeted Next-Generation Sequencing

WGS data was generated for 25 individuals by Complete Genomics (Mountain View, CA, USA) using a mate-paired sequencing technique, as described previously.⁹ Read alignment and variant calling was performed using version 2.5 of the Complete Genomics pipeline.¹⁰

For targeted NGS, enrichments were performed on DNA extracted from peripheral blood using Agilent SureSelect Custom Design target-enrichment kits (Agilent, Santa Clara, CA, USA). Enrichment kits were designed to capture known

pathogenic intronic variants and the protein-coding regions +/-50 nucleotides of selected NCBI RefSeq transcripts; conditions tested included IRD (105 genes or 180 genes), congenital cataracts (114 genes), cardiac disorders (72 genes comprised of 10 sub-panels) and metabolic disorders (226 genes comprised of 6 sub-panels). The genes and transcripts included in the targeted capture regions for each disease referral are available online (Supp Tables S1-S4) and through the UK Genetic Testing Network (<http://ukgtn.nhs.uk/find-a-test/search-by-laboratory/laboratory/manchester-rgc-36/>). Samples were pooled and paired-end NGS was performed using the manufacturer protocols for the Illumina HiSeq 2000/2500 platform (Illumina, Inc., San Diego, CA, USA). Sequencing reads were demultiplexed with CASAVA v.1.8.2. and aligned to the *hg19* reference genome using Burrows-Wheeler Aligner short read (BWA-short v0.6.2) software¹¹ before duplicate reads were removed using samtools v0.1.18. 10.3 million unique NGS reads were generated, on average, per sample ($n=170$, min=1,241,785, max=23,240,481, median=10,812,279), with an average coverage of 880 unique reads per nucleotide ($n=70,514,012$, min=0, max=7956, median=783, sd=515.4) and 2155 unique reads per exon ($n=388,974$, min=0, max=317678, median=1561, sd=3309.8) within the complete region enriched for analysis. The detection and clinical analysis of single nucleotide variants and small insertions/deletions was performed as described previously.^{4, 12}

Copy Number Variant Detection

For the 25 samples with WGS data, CNVs were identified using version 2.5 of the Complete Genomics pipeline.¹⁰ Briefly, this strategy incorporates an assessment of (i) sequencing read depth, and (ii) discordant mate-pairs. For each tested

individual, sequencing read depth was normalized for GC content and genomic positional effects, and CNV status was calculated for non-overlapping 2Kb genomic intervals through comparison to a baseline sample set – comprised of 52 unrelated individuals. To identify the location of breakpoints and insertion points of CNV events, genomic regions where mate-pairs aligned to the reference genome displayed abnormal genomic intervals between the two reads were flagged. Within these flagged regions, local *de-novo* assembly was then performed for sequencing reads where only one of the two reads within the mate-pair aligned to the reference genome. Where possible, the genomic location of breakpoints and insertion points was identified and reported.

For targeted NGS samples, CNV detection was performed using ExomeDepth v1.1.6.¹³ For each tested individual, the ExomeDepth algorithm builds the most suitable reference set from the BAM files of a presented group of potential reference samples. We presented ExomeDepth with BAM files for >20 individuals that had been generated by identical laboratory and computational procedures. All potential reference samples were individuals referred for genomic diagnostic testing who were not knowingly related to the tested individual, and had been obtained from the same sequencing run on the Illumina HiSeq platform, where possible. The reference sample sets selected by ExomeDepth are referred to as ‘reference samples’ herein.

Accuracy of ExomeDepth in comparison to WGS and MLPA

For 25 individuals with IRD we generated gene panel NGS and WGS datasets (Figure 1). We used the variant detection techniques applied to the WGS datasets as a reference standard for CNV detection and then assessed, at the exon level

($n=1590$ exons per sample), the sensitivity and specificity of ExomeDepth applied to gene panel NGS datasets. We defined sensitivity as the capability of ExomeDepth to identify exons with abnormal CNV, and specificity as the capability to identify exons with a normal CNV status. For a further five individuals with cardiac disorders we generated gene panel NGS and MLPA datasets (Figure 1) and then assessed the sensitivity of ExomeDepth applied to gene panel NGS datasets for these individuals.

Assessment of ExomeDepth to identify simulated CNV events

Simulated CNV events were introduced into targeted NGS data for the 25 IRD patients with complementary WGS data. The enrichment region for targeted NGS for the 25 samples encapsulates 1590 protein-coding exons for 105 genes associated with IRD. Importantly, we had previously defined and reported the copy number status for each exon included within the targeted enrichment through the analysis of WGS data.⁷

Simulation was performed using a random sample and exon selector, bedtools v2.16.2 intersect, and software within the PicardTools v1.75 java package: DownsampleSam and MergeSamFiles (Supp Figure S1). Exons were excluded from analysis if they overlapped with known heterozygous deletion events in the selected sample. We simulated deletion events for 1000 single exons and 1749 multiple exons (2, 3 and 4 exon events). In all cases, we assume that the intronic breakpoints of the deletion event are not captured through NGS. Deletion events are not expected to be detected above a test:reference sample read ratio of 0.7 (see supplemental results and methods). We created three discrete groups for

simulated deletion events, with the extent of sequencing reads randomly removed indicated in parentheses: (i) control events (0%), (ii) deletion with amplification bias (40%) and (iii) deletion without amplification bias (50%). Further details on the simulation methodology are provided in the supplemental results and methods.

Assessment of factors influencing successful identification of CNV events

We assessed a number of criteria for known and simulated CNV events in order to assess whether they are key determinants of successful CNV surveillance through ExomeDepth, including: (i) the intra-sample variation in coverage, using the normalized read count (reads-per-kilobase-per-million, *rpk*) coefficient of variation (CV) for surveyed genes in test samples (Supp Figure S2), (ii) the inter-sample variation in coverage, using the *rpk*CV for surveyed exons across reference samples selected by ExomeDepth (Supp Figure S2), (iii) the percentage of nucleotides and the number of exons containing nucleotides with appropriate sequencing depth for in-house diagnostic surveillance (>50x unique sequencing reads), (iv) the total and normalized read depth across surveyed exons, (v) the GC content of the surveyed regions, (vi) the size of exons, and (vii) the distance between neighbouring exons. All statistical analyses were performed in R v3.2.1 software.

Integration of CNV detection during clinical NGS testing

We integrated CNV detection using ExomeDepth into the NGS workflow for 140 individuals from 12 distinct referral groups (Supp Table S5). The reasons for assessment of CNV events were (i) an assessment of whether a heterozygous CNV event was *in-trans* to a clearly or likely pathogenic variant, or (ii) an assessment

of whether a heterozygous CNV event was present in a gene highly specific to an individual's clinical presentation. In accordance with the recommendations of the ExomeDepth developers, test samples with an overall correlation to selected reference samples <0.97 were repeated with an alternative set of reference samples or excluded from analysis. Clinical interpretation of CNVs was restricted to genes relevant to their referral on a case-by-case basis. We performed additional assays to confirm the presence of all identified CNVs before they were clinically reported. Where kits designed and created by MRC-Holland (Amsterdam, Netherlands) were available, we carried out multiplex ligation-dependent probe amplification (MLPA) assays. In the absence of a suitable MLPA kit, we validated CNVs using droplet digital PCR or a bespoke multiplex quantitative fluorescence methodology (see Supplemental Methods). Validated CNV events were submitted to the ClinVar database.¹⁴

Results

Accuracy of ExomeDepth in comparison to WGS and MLPA

To establish the accuracy and reliability of ExomeDepth when applied to targeted NGS data, we analysed targeted NGS datasets for 30 individuals in whom CNV detection had been performed using either WGS ($n=25$) or MLPA ($n=5$). This allowed calculation of the sensitivity and specificity for identified deletions and duplications. Overall, we found a sensitivity of 92.9% and identified that variable and insufficient coverage within surveyed genes reduces the capability of ExomeDepth to identify single exon deletions.

In comparison to WGS, we determined that ExomeDepth applied to targeted NGS datasets (encompassing 1590 exons from 105 genes) has a sensitivity of 100% and a specificity of 99.8% (Supp Table S6) at the exon level. True positive events included a single exon deletion in *GPR98*, a 2 exon deletion in *USH2A*, and a 6 exon deletion in *PCDH15* (Supp Table S7). In comparison to MLPA, we identified 3 out of 4 single exon deletions and one single exon duplication (Supp Table S7). We assessed a number of key factors, and observed that the sequencing data for the individual in whom a single exon deletion was erroneously not identified, showed the highest intra-sample variation (62%) and the highest level of insufficient coverage (9.5% of exons and 0.86% of nucleotides; sample 14011718, Supp Table S8).

We assessed metrics calculated by ExomeDepth for the 8 previously identified deletions and duplication events, observing that the average confidence (Bayes factor, *BF*) determined by ExomeDepth for true positive CNV events was 45.04 (Supp. Table S7, min=6.4, max=76.8) and the average ratio of sequencing reads

between test and reference samples for deletions was 0.61 (Supp. Table S7, min=0.539, max=0.745) and 1.4 for the sole duplication event.

Capability of ExomeDepth to identify simulated CNV events

In order to assess factors that influence the successful identification of CNV events in targeted NGS data using ExomeDepth, we introduced simulated events, *in-silico*, into the targeted NGS datasets created in a clinical setting for the 25 individuals for whom we held complementary WGS data. We found a 97% sensitivity for simulated events when 50% of the NGS reads were removed from selected exons ($n=2749$), and identified that inter-sample variation – a measure of consistency of NGS read coverage across reference samples (Supp Figure S2) – and insufficient coverage were key determinants of whether simulated events were missed or identified by ExomeDepth (Tables 1 & S9).

Single exon deletions ($n=1000$) were introduced into 101 of 105 genes enriched during NGS and we observed that the sensitivity of ExomeDepth for simulated events was 93.5%, with 930 deletions precisely detected at the exon level and 5 included in deletion events erroneously identified as spanning to adjacent exons. This sensitivity is reduced to 79.5% when accounting for amplification bias in simulated events (Supp Tables S10 & S11), with an additional 140 false negative events identified when only 40% of the original NGS reads were removed from the selected exon. Interestingly, 51% (36/70) of the false negative simulated events without amplification bias (50% of NGS reads removed) were exons flanked by neighbouring exons within 250 nucleotides of the canonical donor or acceptor sites. Further, all of these 36 events could be identified if the neighbouring exon boundaries were merged into a single analysis region for simulations, increasing

the overall sensitivity of ExomeDepth for simulated events to 97.1% (Supp Table S11).

Multiple exon deletions ($n=1749$) – where 50% of the NGS reads were randomly removed from adjacent exons – were introduced into all of the 105 genes enriched during targeted NGS for all 25 individuals. We observed sensitivity rates of 96.6% ($n=620$), 95.9% ($n=586$) and 97.1% ($n=543$) for 2 exon, 3 exon and 4 exon deletions, respectively.

To ensure that the process of introducing simulated events into targeted NGS data did not influence the performance of ExomeDepth, we performed the same computational processes of the simulation technique for each event, without removing any NGS reads. No single exon or multiple exon simulated deletion events were identified by ExomeDepth in any of these control simulation experiments.

Integration of CNV detection during clinical NGS testing

Following assessment of the accuracy and the reliability of ExomeDepth applied to targeted NGS datasets, we then integrated CNV detection using ExomeDepth into the NGS workflow for 140 individuals from 12 distinct referral groups to assess specific clinical evaluations. These included either (i) an assessment of whether a heterozygous CNV event was *in-trans* to a clearly or likely pathogenic variant, or (ii) an assessment of whether a heterozygous CNV event was present in a gene highly specific to an individual's clinical presentation. This analysis strategy led to the surveillance of a single gene for 128 individuals, two genes for 10 individuals and three genes for 2 individuals.

Confirmation of molecular diagnoses for 17 individuals

Analysis on a gene-by-patient basis identified 17 heterozygous CNV events (15 deletions, 1 duplication and 1 complex event; Supp Table S12; Supp Figure S3). All events were verified through an alternative technique, were concluded to contribute to the molecular diagnosis for referred individuals and have been submitted to the ClinVar database (Submission number: SUB2171211). The heterozygous CNV events identified by ExomeDepth ranged from a 20 exon deletion in *PCDH15* (NG_009191.2, NM_001142770.1; >600Kb) to single exon deletions in *RPGRIP1* (NG_008933.1, NM_020366.3), *BEST1* (NG_009033.1, NM_004183.3) and *NMNAT1* (NG_032954.1, NM_022787.3). For a single individual referred with a provisional clinical diagnosis of Marfan syndrome, we identified a complex event in *FBN1* (NG_008805.2, NM_000138.4): a 3-exon deletion (chr15:48737523-48741140, c.(5545+1_5546-1)_(5917+1_5918-1)del) and a 2-exon duplication (chr15:48720493-48723049, c.(6739+1_6740-1)_(6997+1_6998-1)dup), consistent with a clinical diagnosis of Marfan syndrome (Figure 2).

We assessed metrics calculated by ExomeDepth for identified deletion and duplication events, observing that the average confidence score (BF) attributed to identified CNV events by the ExomeDepth algorithm was 87 (Supp Table S12, min=22, max=321) and the average read count ratio between test and selected reference samples was 0.56 (min=0.518, max=0.637) and 1.35 (min=1.31, max=1.38), respectively.

Accuracy of ExomeDepth applied in a clinical context

To estimate the accuracy of ExomeDepth applied to targeted NGS datasets for the 123 individuals determined to be absent of CNV events, we assessed (i) copy number variant status through orthogonal techniques, and (ii) two key factors identified through assessments of simulated and known CNV variants: inter-sample variation and insufficient coverage (Table 1).

We calculated the sequencing coverage for each individual, and identified that 3% (135/4551) of the surveyed exons contained at least one nucleotide with less than 50 unique NGS reads. Nine of these exons were found in individuals with a confirmed CNV event in the gene, and 28 were in a gene confirmed to be absent of a CNV event through orthogonal techniques (MLPA; Supp Figure S4). Of the remaining 97 exons, 34 were unique patient-exon combinations and 63 were accounted for by 12 exons with insufficient coverage across multiple samples. On average, 4.6% of the nucleotides within these 97 poor coverage exons received less than 50 unique NGS reads ($n=97$, $\text{min}=0.1\%$, $\text{max}=40.9\%$, $\text{median}=3.6\%$), and all exons were within the range of insufficient coverage values observed for true positive simulated deletion events (Table 1).

To estimate the accuracy of ExomeDepth in relation to reference samples, we calculated the variability of sequencing coverage across the selected references for each individual, and identified an average inter-sample variation for surveyed exons of 5.1% ($n=4551$, $sd=3.4\%$), with average minimum and maximum values observed per-individual of 2.4% ($sd=1.9\%$) and 9.9% ($sd=5.5\%$), respectively. In comparison to simulated single exon deletions, these data are consistent with an average sensitivity of 98.7% ($sd=1.5\%$, $\text{min}=88.7\%$, $\text{max}=100\%$; Figure 3).

For 6 individuals, data from MLPA analyses provided additional support for the absence of a CNV event (Supp Figure S4). For a single individual, we identified a false negative event after subsequent MLPA analysis of the *DSP* gene. We found that alteration of the analysis region, to survey 5 sub-exonic regions enriched by non-overlapping probes though ExomeDepth identified a partial exon duplication event within the *DSP* gene which complemented the result from MLPA (Supp Figure S5).

Discussion

Copy number variants (CNVs) are an important and common form of genomic variation in the general population,^{15, 16} and are implicated in many Mendelian disorders.^{7, 8, 17} An ability to accurately survey for CNV events, in particular in targeted NGS datasets, therefore has the power to increase diagnostic yields and enhance clinical care. While it has already been shown that read count CNV detection algorithms can be successfully applied to targeted NGS data in a research context,^{13, 18-20} their integration within diagnostic services has been slower due to a lack of validation parameters. In this study, we have identified key factors which can facilitate the successful application of a widely used bioinformatics tool, ExomeDepth,¹³ for CNV surveillance of targeted NGS datasets within the clinical environment.

CNV detection tools used in a diagnostic context must be able to identify deletion and duplication events that encapsulate single targets/exons included within the targeted enrichments of custom gene panel and WES techniques, which is a known limitation of some publically available algorithms. Since large datasets of known true positive single exon CNV events do not exist, we have developed and applied a computational simulation technique which permits extended assessment of single exon CNV events. As a result, we have been able to perform an assessment of trends in large and controlled datasets (Table 1), We have then used real-time comparison between WGS and targeted NGS data to assess their applicability to real datasets. Using this combined approach we have shown that amplification bias within NGS assays and the distance between exons enriched during NGS influences the overall sensitivity of ExomeDepth (Supp Table S11). After accounting for these dominating factors, we have demonstrated how variability of sequencing coverage between and within samples, the extent of read

depth, the size of surveyed exons and the level of insufficient coverage are important determinants of successful identification of single exon deletion events through ExomeDepth (Table 1, Table S9 and S10). Whilst all these metrics are indicated as important quality assurance parameters for the accurate detection of single exon CNVs, they are neither completely independent nor equally applicable to real datasets on an individual basis. We therefore selected two key metrics for routine incorporation into diagnostics: insufficient coverage (test sample dependent) and inter-sample variability (reference sample dependent). This two-part process firstly checks for the quantity of sequencing coverage over exons surveyed in the tested sample, and second, assesses the consistency of NGS read coverage across reference samples for each surveyed exon. We have assimilated this information to successfully integrate surveillance of CNVs into the clinical bioinformatics pipeline for 140 individuals in a clinical setting, achieving a definitive molecular diagnosis in 17 of 140 individuals. Importantly, we have shown that 97.2% of the exons surveyed and determined to be absent of a CNV event have sufficient coverage, and none of the insufficiently covered exons lie outside the range of true positives identified from simulated experiments. Moreover, we have calculated the inter-sample variability for surveyed exons on an individual basis, and through comparison to simulated single exon events, estimated the accuracy of ExomeDepth to be 98.7% for the 123 individuals without an identified CNV (Figure 3). Both of these quality assurance observations are supported by their integration with other CNV software tools²¹ and the absence of CNV events in 6 individuals tested through MLPA.

Taken together, our data illustrate the utility of CNV assessments within a diagnostic setting using the publically available ExomeDepth software, and support the utilization

of quality assurance parameters in complement to CNV detection algorithms in targeted NGS diagnostic services. Whilst other types of software can be routinely applied to WGS datasets to detect CNVs at single nucleotide resolution, we expect that application of the approaches outlined in this study will improve the utilization of read depth CNV tools in diagnostic environments across heterogeneous targeted NGS gene panel approaches, including small and large gene panels, as described here, and WES.

Declarations

Competing interests

The authors of this manuscript have no competing interests to declare.

Funding

This work was supported by the Biotechnology and Biological Sciences Research Council, the Manchester Biomedical Research Centre, The National Institute for Health Research Biomedical Centre at Moorfields Eye Hospital and the UCL Institute of Ophthalmology, the DST UK-India Education and Research Initiative, and an independent research grant funded by the Manchester Academic Health Science Centre.

Acknowledgments

We thank all patients, referring clinicians, clinical scientists and genetic counselors involved in this study.

Authors' contributions

JME, GCMB, SCR and BP designed and coordinated the study. JME, CC, StB, SaB, SG, RLG, PIS, BH, MM, ARW, WGN, SCR, GCMB contributed genetic and/or phenotypic data. JME wrote the manuscript and all authors provided important revisions and intellectual content.

Bibliography

1. Baetens M, Van Laer L, De Leeneer K, et al. Applying Massive Parallel Sequencing to Molecular Diagnosis of Marfan and Loeys-Dietz Syndromes. *Human Mutation* 2011; **32**: 1053-62.
2. O'Sullivan J, Mullaney BG, Bhaskar SS, et al. A paradigm shift in the delivery of services for diagnosis of inherited retinal disease. *J Med Genet* 2012; **49**: 322-6.
3. Nishio SY, Hayashi Y, Watanabe M, Usami SI. Clinical Application of a Custom AmpliSeq Library and Ion Torrent PGM Sequencing to Comprehensive Mutation Screening for Deafness Genes. *Genet Test Mol Biomarkers* 2015.
4. Ellingford JM, Barton S, Bhaskar S, et al. Molecular findings from 537 individuals with inherited retinal disease. *J Med Genet* 2016.
5. Lee H, Deignan JL, Dorrani N, et al. Clinical Exome Sequencing for Genetic Identification of Rare Mendelian Disorders. *JAMA* 2014.
6. Yang Y, Muzny DM, Xia F, et al. Molecular Findings Among Patients Referred for Clinical Whole-Exome Sequencing. *JAMA* 2014.
7. Ellingford JM, Barton S, Bhaskar S, et al. Whole Genome Sequencing Increases Molecular Diagnostic Yield Compared with Current Diagnostic Testing for Inherited Retinal Disease. *Ophthalmology* 2016; **123**: 1143-50.
8. Gilissen C, Hehir-Kwa JY, Thung DT, et al. Genome sequencing identifies major causes of severe intellectual disability. *Nature* 2014.
9. Drmanac R, Sparks AB, Callow MJ, et al. Human genome sequencing using unchained base reads on self-assembling DNA nanoarrays. *Science* 2010; **327**: 78-81.

10. Carnevali P, Baccash J, Halpern AL, et al. Computational techniques for human genome resequencing using mated gapped reads. *J Comput Biol* 2012; **19**: 279-92.
11. Li H, Durbin R. Fast and accurate short read alignment with Burrows-Wheeler transform. *Bioinformatics* 2009; **25**: 1754-60.
12. Gillespie RL, O'Sullivan J, Ashworth J, et al. Personalized diagnosis and management of congenital cataract by next-generation sequencing. *Ophthalmology* 2014; **121**: 2124-37.
13. Plagnol V, Curtis J, Epstein M, et al. A robust model for read count data in exome sequencing experiments and implications for copy number variant calling. *Bioinformatics* 2012; **28**: 2747-54.
14. Landrum MJ, Lee JM, Benson M, et al. ClinVar: public archive of interpretations of clinically relevant variants. *Nucleic Acids Res* 2016; **44**: D862-8.
15. Iafrate AJ, Feuk L, Rivera MN, et al. Detection of large-scale variation in the human genome. *Nat Genet* 2004; **36**: 949-51.
16. Sebat J, Lakshmi B, Troge J, et al. Large-scale copy number polymorphism in the human genome. *Science* 2004; **305**: 525-8.
17. Glessner JT, Bick AG, Ito K, et al. Increased frequency of de novo copy number variants in congenital heart disease by integrative analysis of single nucleotide polymorphism array and exome sequence data. *Circ Res* 2014; **115**: 884-96.
18. Krumm N, Sudmant PH, Ko A, et al. Copy number variation detection and genotyping from exome sequence data. *Genome Res* 2012; **22**: 1525-32.

19. Fromer M, Moran JL, Chambert K, et al. Discovery and statistical genotyping of copy-number variation from whole-exome sequencing depth. *Am J Hum Genet* 2012; **91**: 597-607.
20. Li J, Lupat R, Amarasinghe KC, et al. CONTRA: copy number analysis for targeted resequencing. *Bioinformatics* 2012; **28**: 1307-13.
21. Johansson LF, van Dijk F, de Boer EN, et al. CoNVaDING: Single Exon Variation Detection in Targeted NGS Data. *Hum Mutat* 2016; **37**: 457-64.

Figure Legends

Figure 1. Study Design. The approach taken in this study to assess the 'accuracy' and 'key factors' influencing the accuracy of ExomeDepth applied to targeted next-generation sequencing datasets. The 'key factors' assessed by application of ExomeDepth to datasets with known and simulated CNVs are outlined in Table 1. *CNV*, copy number variation; *WGS*, whole genome sequencing; *MLPA*, multiplex ligation-dependent probe amplification; *gene panel*, next-generation sequencing data generated in a diagnostic environment after enrichment for a set of genes known as a cause of specific Mendelian disorders.

Figure 2. *FBN1* copy number variant. A complex 3-exon deletion, $c.(5545+1_5546-1)_{(5917+1_5918-1)}$ del, and 2-exon duplication, $c.(6739+1_6740-1)_{(6997+1_6998-1)}$ dup event identified in *FBN1* (NG_008805.2, NM_000138.4), confirming a clinical diagnosis of Marfan syndrome for the referred individual. *Red crosshairs*, the ratio of reads between test and reference samples; *grey bar*, the 95% confidence interval of expected read ratios in comparison to reference samples.

Figure 3. Inter-sample variation in sequencing coverage across surveyed exons. *Simulations*, the variability of sequencing coverage in selected reference samples for 971 identified and 29 missed single exon simulated deletions. *Diagnostic survey*, the variability of sequencing coverage in selected reference samples for 4551 exons surveyed for copy number variants in a diagnostic context.